| | | |
|---------------|----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| INERIS | PROTOCOLE DE MESURE DES COMPOSES [A], [B], [C], ETC. DANS LES PAPIER THERMIQUES | MEASUREMENT PROTOCOL OF COMPOUNDS [A], [B], [C], ETC. IN THERMAL PAPERS |
|---------------|----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|

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This information document describes the measurement protocol of compounds [A], [B], [C]... in thermal papers. Each appendix refers to the operating method of the particular compound(s) to be measured. It also gives, on an exploratory basis, the threshold(s) below which thermal paper is considered exempt of the concerned compounds [A], [B], [C].

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LIST OF APPENDICES

Appendix 1: operating method for the analysis of BISPHEOLS A, F and S in thermal papers

Appendix 2: operating method for the analysis of 15 phenolic compounds, including BISPHEOLS A, F and S, in thermal papers

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1. Scope:

This operating method describes the dosage of bisphenols A, F and S in thermal papers by liquid chromatography. It describes:

- a single extraction method;
- three analytical methods that can be used, at the laboratory's discretion: liquid chromatography with UV detection (diode arrays), liquid chromatography with fluorescence detection, liquid chromatography couple to a tandem mass spectrometer.

Whatever the analytical method used, limits of quantification lower or equal to 0.01% (100 µg/g of paper) can be reached for each bisphenol.

2. Principle:

Bisphenols are extracted from thermal paper by methanol at ambient temperature. The three bisphenols are separated using liquid chromatography on a C18 column. They are detected and quantified by UV detection (diode arrays), fluorescence, or mass spectrometry.

3. Precautions:

Thermal papers must be conserved away from light and without direct contact between samples (for example in a glass, PP or PE bottle, in aluminum foil...) to avoid any cross contamination: migration of bisphenols between samples is observed, especially when thermal

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papers are slightly wet.

Whatever is the detection method used, the observed interferences for thermal paper analysis are limited. Nevertheless, detection by mass spectrometry allows to reduce the risk of interference.

4. Extraction:

The extraction method is performed at ambient temperature:

- 50 mg of paper are introduced into a glass, PE or PP vial;
- 5 ml of methanol are added into the vial: the paper must be completely immersed in methanol;
- The vial is placed under vortex agitation for one minute;
- The vial is then placed for 10 minutes in an ultrasonic bath.

Extracts must be stored at 4°C in a dark refrigerator before analysis.

Note: For thermal papers thicker than receipts, a test of at least 400 mg of the sample is advised. The paper must first be cut into small pieces. The volume of methanol is increased (10 ml) in order to ensure that the thermal paper is completely immersed in the solution.

5. Quantitative Determination:

The bisphenols are analysed by liquid chromatography. They can be detected and quantified by a diode array detector, a fluorimeter or a mass spectrometer. These three techniques allow to reach limits of quantification less or equal to

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0.01% (100 µg/g of paper) and they are not subject to notable interferences during the application of this procedure. Each laboratory can thus choose the analytical technique that seems the most appropriate in relation to the disposal material and its expertise.

5.1 Liquid Chromatography Coupled with a Diode Array Detector:

The chromatographic conditions are the following:

- C18 column 250mm*4.6mm 5µm
- Eluent: methanol/10mM of H₃PO₄ in water (60/40)
- Flow rate of 1ml/min
- Isocratic conditions for 20 min then 10 min at 100% methanol to clean the column
- Injection volume of 20µl
- Diode array detector: 230 nm
- Quantification by external calibration from 0.5 to 200µg/ml

The use of a simple UV detector is not enough to confirm the compounds identity. It's preferred to work with a diode array detector to confirm components by comparing their absorption spectrum. The following quantification limits were validated for this technic:

- Bisphenol A: 0.005% (50 µg/g of paper),
- Bisphenol F: 0.005% (50 µg/g of paper),
- Bisphenol S: 0.005% (50 µg/g of paper)

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Figure 1: chromatogram obtained by HPLC - diode arrays (230 nm)

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5.2 Liquid Chromatography Coupled to a Fluorimeter:

The chromatographic conditions are the following:

- C18 column 250mm*4.6mm 5µm
- Eluent: methanol/10mM of H₃PO₄ in water (60/40)
- Flow rate of 1ml/min
- Isocratic conditions for 20 min then 10 min at 100% methanol to clean the column
- Injection volume of 20µl
- Excitation wave length: 275 nm
- Emission wave length: 313 nm
- Quantification by external calibration from 0.01 to 1µg/ml for BPF and BPA, from 1 to 100µg/ml for BPS

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Bisphenols A and F are very fluorescent:

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the quantification limits that can be reached are better than that measured in UV detection. However, bisphenol S is not very sensitive for fluorescence and the obtained LOQ is comparable to that determined by UV. The following quantification limits were validated:

- Bisphenol A: 0.0001% (1 µg/g of paper),
- Bisphenol F: 0.0001% (1 µg/g of paper),
- Bisphenol S: 0.01% (100 µg/g of paper)

5.3 Liquid Chromatography Coupled to a Tandem Mass Spectrometer

The chromatographic conditions are the following:

- C18 X-Bridge Column
2.1mm*50mm 2.5µm
- Eluent A: 2 mM of ammonium acetate in methanol
- Eluent B: 2 mM of ammonium acetate in water
- Flow rate: 0.3ml/min
- Gradient:

| Time in minutes | %A | %B |
|-----------------|-----|----|
| 0 | 10 | 90 |
| 10 | 65 | 35 |
| 11 | 100 | 0 |
| 12 | 100 | 0 |
| 13 | 10 | 90 |
| 15 | 10 | 90 |

- Injection volume is 3µl
- Operating parameters of the TQD (Triple Quadropole from Waters®):
 - Polarity: ES-

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-Quantification by external calibration from 0.5 to 200 µg/ml

Limits of quantification equivalent to those obtained in UV and fluorescence are validated. Lower quantification limits can be easily reached with this analytical technique but this

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optimization is not necessary for this application.

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Figure 2: example of a chromatogram obtained by LC-MS

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Note: Other types of mass spectrometers can be used if the operatory conditions have been validated.

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6. Quality Controls:

Quality control is guaranteed by the systematic performance of a method blank, method control and drift control.

The **method blank** (extraction without paper) must be performed to guarantee the absence of contaminants coming from the reagents or vials used during the analysis: the content of each phenolic compound in the paper must be under 0.01%.

The **method control** allows to verify that extraction yield for each component is satisfactory (recovery rates greater than 90%). This control is performed by adding a calibration solution on the thermal paper before extraction (50 µl of solution on 50 mg

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of paper).

Drift control must be performed throughout the analytical sequence to guarantee the absence of drift during the analytical sequence. This monitoring is performed by injecting a calibration solution at a regular frequency (for example every ten samples). This drift must be limited to 15% during the whole analytic sequence.

7. Content thresholds:

It is given, on an exploratory basis, that the thermal paper is considered to be exempt of bisphenol A when the content threshold is less than 0.02% by mass ($2 \cdot 10^{-4}$ by mass).

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1. Scope

This operating method describes analysis techniques for 15 phenolic compounds, including bisphenols A, F and S, in thermal papers. These phenols are used as additive on the thermo-sensitive coating for thermal paper. The 15 molecules subject of this operating method are (see chemical formulas at the end of the document):

- Bisphenol A (BPA) CAS 80-05-7
- Bisphenol S (BPS) CAS 80-09-1
- Bisphenol F (BPF) CAS 620-92-8
- 2,4'-Bisphenol S (2,4'-BPS) CAS 5397-34-2
- TGSA CAS 41481-66-7
- DD70 CAS 93589-69-6
- D8 CAS 95235-30-6
- Bisphenol AP (BPAP) CAS 1571-75-1
- Bisphenol C (BPC) CAS 79-97-0
- Bisphenol PH (BisOPP-A) CAS 24038-68-4
- BPS-MAE CAS 97042-18-7
- BPS-MPE CAS 63134-33-8
- ISO F CAS 2467-02-9
- PHBB CAS 94-18-8
- D90 CAS 191680-83-8

This operating method describes:

- a single extraction method;
- two analytical methods that can be used at the laboratory's discretion:
 - liquid chromatography with UV

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détection (diode arrays) and confirmation by fluorescence detection,

- liquid chromatography coupled to tandem mass spectrometry.

Whatever the analytical method used, limits of quantification lower or equal to 0.01% (100 µg/g of paper) can be reached for each bisphenol.

2. Principle

Phenolic compounds are extracted from thermal paper by methanol at ambient temperature. The three bisphenols are separated using liquid chromatography on a C18 column. They are detected and quantified either by UV detection (diode arrays) coupled to a fluorescence detector, or by tandem mass spectrometry.

3. Precautions

Thermal papers must be conserved away from light and without direct contact between samples (for example in a glass, PP or PE bottle, in aluminium foil...) to avoid any cross contamination: migration of bisphenols between samples is observed, especially when thermal papers are slightly wet. Whatever is the detection method used, the observed interferences for thermal paper analysis are limited. Nevertheless, detection by mass spectrometry allows to reduce the risk of interference.

4. Extraction:

The extraction method is performed at ambient temperature:

- 50 mg of paper are introduced into a glass, PE or PP vial;

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- 5 ml of methanol are added into the vial: the paper must be completely immersed in methanol;

- The vial is placed under vortex agitation for one minute;

- The vial is then placed for 10 minutes in an ultrasonic bath.

Extracts must be stored at 4°C in a dark refrigerator before analysis.

Note: *For thermal papers thicker than receipts, a test of at least 400 mg of the sample is advised. The paper must first be cut into small pieces. The volume of methanol is increased (10 ml) in order to ensure that the thermal paper is completely immersed in the solution.*

5. Quantitative determination:

The bisphenols are analyzed by liquid chromatography. Two types of detection can be used: either UV detection with a diode array detector (confirmation with a fluorimeter is recommended) or by mass spectrometry. These two techniques allow to reach limits of quantification less or equal to 0.01% (100 µg/g of paper) and they are not subject to notable interferences during the application of this procedure. Each laboratory can thus choose the analytical technique that seems the most appropriate in relation to the disposal material and its expertise.

a. Liquid chromatography coupled with a diode array detector:

Separation is performed on an UHPLC with the following chromatographic conditions:

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- Eluent A: acétonitrile
- Eluent B: 10 mM of H3PO4 in water
- Flow rate: 0.6 ml/min
- Isocratic conditions 12 minutes until 47% acétonitrile, then 12 minutes until 100% acétonitrile and finish with two minutes at 100% to clean the column before returning to initial conditions.

| time (minutes) | Eluant A | Eluant B |
|----------------|----------|----------|
| 0 | 35 | 65 |
| 12 | 47 | 53 |
| 24 | 100 | 0 |
| 26 | 100 | 0 |
| 26.5 | 35 | 65 |
| 30 | 35 | 65 |

- Injection volume: 10 µl
- Diode array detector: 230 nm
- Quantification by external calibration from 0.5 to 200 µg/ml in methanol.

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Figure 1: chromatogram obtained by HPLC - diode arrays (230 nm) for phenolic compounds (other than D90)

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The use of a UV diode array detector allows to verify the presence of compounds by comparing their spectrum to those of standard compounds. To ensure the selectivity of analyses and avoid the risk of a false positive, it is recommended to confirm the results by fluorescence detection (295 nm) for the following phenolic compounds: BPF, BPA, ISO F, PHBB, BPC and BPAP. However, fluorescence detection is not suitable for the other compounds, of which the fluorescence response is weak or insignificant.

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A limit of quantification of 0.005% is validated for all compounds, which corresponds to 50 µg/g of paper.

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b. Liquid chromatography coupled to a tandem mass

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spectrometer

Tandem mass spectrometry has greater sensitivity and better analytical selectivity. The use of ionic fragments to identify a compound gives precise information on the identity of the detected compound. The chromatographic conditions are similar to those used for UV detection, except that ammonium acetate is used instead of phosphoric acid.

- X-Bridge Shield RP-18 column, 3.5µm, 3.0*150mm
- Eluent A: acetonitrile
- Eluent B: 2 mM in ammonium acetate in water
- Flow rate: 0.6ml/min
- Injection volume: 3 µl

| Time (minutes) | Eluent A | Eluent B |
|----------------|----------|----------|
| 0 | 35 | 65 |
| 12 | 47 | 53 |
| 24 | 100 | 0 |
| 26 | 100 | 0 |
| 26.5 | 35 | 65 |
| 30 | 35 | 65 |

Each compound is identified by 2 daughter ions, after fragmentation of the parent ion. All the compounds are analysed in ESI - mode:

| Compound | Molar mass | Parent ion (m/z) | Daughter ions (m/z) |
|----------|------------|------------------|---------------------|
| | | | |

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Note: Other types of mass spectrometers can be used on the condition that the operatory conditions have been validated.

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c. Special case of D90

D90 is less soluble in methanol comparing to the other phenolic compounds. For the calibration solution, it is necessary to prepare the stock solution of D90 in DMSO

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before performing dilutions in methanol. The maximum concentration of DMSO in the solutions injected into the chromatographic column must not exceed 2%. Due to its chemical structure, the chromatogram corresponding to D90 (see figure 3) is contain several peaks, corresponding to several isomers (n = 1, 2, ...). The 2,4'-bisphenol S is an impurity contained in D90 standard.

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Figure 2: Chemical structural formula of D90

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Because of his low solubility in methanol, the extraction yield of D90 in thermal paper using methanol is around 80%.

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Figure 3: Chromatogram obtained by HPLC -

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diode arrays (230 nm) for D90

6. Quality controls:

Quality control is guaranteed by the systematic performance of a method blank, method control and drift control.

The **method blank** (extraction without paper) must be performed to guarantee the absence of contaminants coming from the reagents or vials used during the analysis: the content of each phenolic compound in the paper must be under 0.01%. The **method control** allows to verify that extraction yield for each component is satisfactory (recovery rates greater than 90%). This control is performed by adding a calibration solution on the thermal paper before extraction (50 µl of solution on 50 mg of paper).

Drift control must be performed throughout the analytical sequence to guarantee the absence of drift during the analytical sequence.

This monitoring is performed by injecting a calibration solution at a regular frequency (for example every ten samples). This drift must be limited to 15% during the whole analytic sequence.

7. Content thresholds:

It is granted, on an exploratory basis, that the thermal paper under consideration does not have added phenolic compounds when the content threshold of each of these 15 phenolic compounds is less than 0.02% by mass ($2 \cdot 10^{-4}$ by mass).

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